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 **Terms and conditions of use.**

Disclaimer

Carefully read the present instructions for the use of **Tree•Oss®** dental implants. The effectiveness of the use of our products depends to a great extent on the meticulous observance of the guidelines that will be provided by the intervening professional. This product is part of an integral system and should only be used in conjunction with the original associated products according to **Tree•Oss®** instructions and recommendations. The non-recommended use of products manufactured by third parties will void any warranty or other obligation, express or implied. The user is responsible for determining whether or not a product is indicated for each patient and each circumstance. **Tree•Oss®** disclaims any liability, whether implicit or explicit, and will not be liable for any direct, indirect, punitive or other damage that is caused by or related to any error in the trial or the practice of the professional in the use of **Tree•Oss®** products. The use of this product is under the control of the user, so this will be your responsibility. **Tree•Oss®** does not assume any responsibility for the damages produced as a result.

Description and intended use

Tree•Oss® dental implants are manufactured from grade 4 biocompatible titanium and

sterilized with the aim of replacing lost dental pieces, through dental surgery, both in the upper and lower jaw allowing fixed oral rehabilitation, both individual and multiple, as in turn full or partial removable prostheses to act as retentive support.


Sterilization

Tree•Oss® dental implants are marketed in sterile condition as indicated on the exterior labeling. The sterilization is performed by gamma rays.

Do not use **Tree•Oss®** dental implants whose packaging is damaged or broken. Do not remove **Tree•Oss®** dental implants from their packaging until the moment of use. Under no circumstances re-sterilize.

The strict respect of the chain of custody of the product and adequate asepsis by the professional in the clinical environment, are essential for the product to offer its skills in its entirety.

PACKAGING DESCRIPTION

Storage conditions

Tree•Oss® dental implants have a shelf life of five years, from the moment they are manufactured, under normal conditions of storage in cool places and with low humidity. The product must be stored in its original unopened packaging and in a dry place with temperatures between 5 ° C to 50 ° C. Do not use dental implants after the expiration date.

Precautions and warnings

The implantation of **Tree•Oss®** dental implants must be carried out by professionals who have carried out specific courses in dental implantology in order to avoid the use of inadequate surgical techniques and an incorrect evaluation of the residual bone volume in the jaws, preventing early

failures in the osseointegration as well as bone loss in the prosthetic load. Under no circumstances should **Tree•Oss®** dental implants be used by professionals not trained for this purpose.

If the professional does not respect the limitations of use or the indicated work steps, failures could occur, which can not guarantee the success of the implant at 100%. Implant treatment can result in bone loss or biological or mechanical failure, including fracture of the implants due to fatigue. The close collaboration between the surgeon, the restorative dentist and the dental laboratory technician is essential for the success of the implant treatment.

The osteotomy in the area to be implanted should be done following the recommended

sequence in the standard surgical protocols. It is strongly recommended the constant irrigation of the strawberry during surgical procedures. It is suggested to always use **Tree•Oss®** dental implants, instruments and prosthetic elements designed by **Tree•Oss®**.

It must be taken into account, when removing the **Tree•Oss®** dental implant from its container, that it is not convenient to contact the surface of the implant with alloys other than titanium, to avoid alteration and / or contamination of it. It is also recommended that the **Tree•Oss®** dental implant be placed as quickly as possible once removed from its container to prevent container to prevent adhesion of particles normally found in the air, which can have an impact on the correct osseointegration and create risks of infection, in both in this way the sterilization of the product is affected.

Special attention should be paid to patients presenting local or systemic factors that may affect the healing process, either bone or soft tissue, or the process of osseointegration (for example, if the patient is a smoker, does not maintain proper oral hygiene, suffers from uncontrolled diabetes, is undergoing orofacial radiotherapy or steroid therapy, or has infections in the adjacent bone).

In general, the placement of the implant and the prosthetic design should be adjusted to the individual conditions of the patient. In case of bruxism or unfavorable intermaxillary relationships, a reassessment of the treatment option should be considered. A deficit of hard or soft tissue before the intervention can cause an undesirable esthetic result or unfavorable implant angulations.

In case the **Tree•Oss®** dental implant is placed using the indicated technique for two surgeries (submerged implant), before making the last suture, it must be certain that the interior of the implant is clean of blood residues.

Due to the small size of the devices, attention must be paid so that the patient does not swallow or aspirate them.

Tree•Oss® dental implants are not designed to be reused, as it is indicated in the external labeling, because blood and bone remnants from previous uses would generate an immediate rejection by the organism of the new patient.

The professionals must ensure that the patient is correctly informed about the procedure and the stages corresponding to the postoperative period. It is necessary that the patient is correctly informed, if possible by written means, about the necessary care and the medication to be provided after the surgical procedure.

Never exceed an insertion torque of 45 Ncm during implant placement. Excessive tightening of the implant can cause damage to the implant, fracture or necrosis of the bone bed. If a surgical instrument is used to insert the implant, special care should be taken not to overtighten it.

Implants placed in the jaw should not pierce the membrane of the sinus floor.

After implant placement surgery, in order to guarantee a good long-term treatment result, it is recommended to carry out a complete periodic follow-up of the patient, as well as inform him of the appropriate oral hygiene procedures.

The prosthetic rehabilitation on **Tree•Oss®** implants must respect the basic biomechanical parameters checked in the literature. Take into account the appropriate number of implants placed according to the rehabilitated parts, the maximum lever arm, the occlusion and the correct alignment and distribution of the forces. Avoid large lateral loads or bridges on implants, either between implants or in cantilevers. In turn, a treatment plan with an oversized number of implants can also cause the loss of them. Do not use implants with a diameter smaller than 4mm in the posterior region. Narrow diameter implants, such as 3.5mm or 3.7mm, are not recommended for molar rehabilitation.

Contraindications

A planning and evaluation of the clinical case must be carried out before the surgical procedure to avoid the patient any risk associated with medical conditions that may turn into absolute or relative contraindications for surgery. In the first case, surgery should not be performed while in the second case preventive measures should be taken to avoid the failure of **Tree•Oss®** dental implants. **Tree•Oss®** dental implants should not be placed in patients who are not in the conditions necessary to undergo any general oral surgery. Do not use implants when the patient is being treated with bisphosphonate.

Tree•Oss® implants should not be placed nor should they remain in function if the volume of healthy alveolar bone is not sufficient to support the implant minimally (minimum: 1 mm circumference around the entire implant and 2 mm apex) or if the minimum distances are not respected with a tooth or other implant (2 mm minimum distance with a natural tooth and 3 mm minimum distance with another implant)

Poor bone quality, inadequate oral hygiene of the patient, tobacco abuse, poorly controlled systemic diseases (diabetes, etc.), immunodeficiency, alcoholism, drug addiction and psychological instability may contribute to a lack of osseointegration and the subsequent failure of the implant, both immediately after surgery and after achieving initial osseointegration. A severe bruxism, a prosthetic overload, the loss of cervical bone or the excessive inclination of the implants that produce strong forces on the same off its axis, can cause bone loss, screw loosening, rupture or implant failure. Exposure to radiation and chemotherapy can affect the health of the implant. Dental implant patients should be instructed to consult their doctor before undergoing these types of treatments.

Possible contraindications for dental surgery are:

- Alterations in blood coagulation
- Systemic or local oral infections

- Diabetes
- Smoking
- Patients subject to chemotherapy
- Alcoholism
- Kidney conditions
- Diseases that modify the bone metabolism
- Coronary or vascular conditions untreated
- Insufficient residual bone
- Bruxism and related conditions
- Allergic or hypersensitive to commercially pure titanium (grade 4).

For more information on contraindications, please consult clinical manuals and corresponding textbooks.

Adverse reactions

No literature or scientific publications have been registered that report on the rejection of titanium alloys by the human body, or collateral effects derived from the use of titanium. Failures occur during oral rehabilitation due to primary infections, contamination or incorrect sterilization of the instruments used during the surgical procedure, by the placement of the dental implant in areas with insufficient residual bone or by bone losses derived from the inadequate location of the dental implant with respect to the prosthetic rehabilitation performed. The quality of the **Tree•Oss®** implants, together with a correct observation of the intervening professional of the positioning and manipulation techniques that arise from the supra provided guidelines, guarantees the optimal functionality of the **Tree•Oss®** dental implants.

Warranty

The guarantee offered by **Tree•Oss®** on its dental implants can only be made effective in situations of strict compliance by the professionals involved in these conditions of use.

Safety and Performance

The summary of safety and clinical performance (SSCP) is available in the European Database on Medical Devices (EUDAMED). It can be accessed using the Basic UDI-DI:

Product	Basic UDI-DI
Anatomical IH Dental Implant	84357129IDAIHNC
HS IH Dental Implant	84357129IDHSIHQ7

You can access the European Database on Medical Devices (EUDAMED) at the following link:

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