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<https://ec.europa.eu/tools/eudamed>

Abutments and prosthetic accessories - Instructions for Use

Important: read carefully

Disclaimer:

These prosthetic components are part of a global concept and belong to an integrated dental implant system. They should only be used in conjunction with the original associated products according to **Tree•Oss®** instructions and recommendations. The use of products manufactured by third parties together with the prosthetic components **Tree•Oss®** will void any warranty or other obligation, implicit or explicit, by **Tree•Oss®**. The user of prosthetic components **Tree•Oss®** is responsible for determining whether or not a product is indicated for each patient and each particular clinical case. **Tree•Oss®** declines any responsibility, whether implicit or explicit and will not be liable for any direct, indirect, punitive or other damage that occurs or is related to any error in the judgment or practice of the professional in the use of the **Tree•Oss®** products. The professional is also obliged to regularly study the latest developments and applications with regard to **Tree•Oss®** products. In case of doubt, the professional should contact **Tree•Oss®**. The processing and surgical application of this product are under the control of the intervening professional and this will be your responsibility. **Tree•Oss®** does not assume any responsibility derived from the damages produced as a consequence of its incorrect use.

With the instructions for cleaning and sterilization, **Tree•Oss®** provides a validated procedure to guarantee sterile products. It is guaranteed that the conditioning will be carried out using appropriate equipment, materials, and personnel, in suitable facilities to ensure the desired result. It is the user's responsibility. In the same way, it is necessary for the user to correctly evaluate any deviation that is made regarding the instructions supplied, in order to determine the effectiveness and the possible adverse consequences.

Description:

The **Tree•Oss®** prosthesis family is comprised of:

- ◆ Dental abutments
- ◆ Healing caps

Prefabricated prosthesis for endosseous dental implants with a direct connection to the dental implant, designed for use in prosthetic rehabilitation.

Indication:

Tree•Oss®'s prosthesis are indicated for cemented and screwed crowns of unitary teeth, as well as for bridges of multiple units cemented on implants.

Contraindications:

Used in clinical cases indicated correctly, it has no contraindications.

Instructions for clinicians:

We recommend that clinicians, both new users and users with experience in implants, receive specialized training before using a new treatment method. For more information visit www.treeossiberica.es. Always work with an experienced professional the first time you use a new treatment method.

Precautions during the procedure:

The close collaboration between the surgeon, the dentist and the laboratory technician is essential to achieve a successful implant treatment.

It is especially important to achieve a correct distribution of forces by adapting and passively adjusting the bridge to the abutments of the implant by adjusting the occlusion to the opposite jaw and avoiding excessive lateral load forces, especially in cases of immediate loading. Due to the small size of the prosthetic components, care must be taken so that the patient does not swallow or aspirate them.

Storage

Store the **Tree•Oss®** prosthesis in a dry, clean, dust-free environment at moderate temperatures of 5 °C to 50 °C (41 °F to 122 °F). To prevent contamination, keep the prosthesis in their original container.

For information on surgical procedures and protocol, refer to the **Tree•Oss®** implant instructions at www.treeossiberica.es.

Useful life:

Under normal conditions of use and subject to non-harmful occlusal forces, the useful life of the **Tree•Oss®** Abutment is 10 years.

Sterility:

Tree•Oss® prosthesis and Prosthetic Accessories are supplied clean and non-sterile. Before using the product for the first time you must sterilize it. To do this place a single device in a sterilization bag, seal the bag and autoclave by steam, respecting the following parameters:

- ♦ Temperature: 134 °C (273,2 °F)
- ♦ Sterilization time: 5 minutes
- ♦ Drying time: 40 minutes

When loading the autoclave with more than one bag, place them in an upright position, side by side, with the transparent sheet surfaces facing each other. After sterilization, place the devices in a dry, clean and dust-free place. Follow the bag manufacturer's instructions regarding the storage conditions and the expiration date of the sterilized products.

Disposal:

These elements are arranged as pathogenic residues after their contact with the patient.

Caution:

The use of non-sterile prosthesis can cause infectious diseases or infection of tissues. This product can only be sold to dentists or authorized doctors. The product is designed for single use only. An abutment that has been in contact with one patient cannot be reused with another patient, as this carries the risk of cross-contamination.

Precaution:

The use of non-sterile prosthesis may cause infectious diseases or tissue infection. The product is designed for single use only. A prosthesis that has been in contact with one patient must not be reused with another patient, as this poses a risk of cross-contamination.

This product may only be sold to licensed dentists or authorized physicians.

Safety and Performance:

The Summary of Safety and Clinical Performance (SSCP) is available in the European Database on Medical Devices (EUDAMED) and can be accessed through the Basic UDI-DI:

Product	Basic UDI-DI
Straight Standard Abutment Ti	84357129PRETiX8
Straight Abutment without Shoulder Ti	84357129PRSHTi3E
Angled Esthetic Abutment Ti	84357129PEsAng5L
Multiple abutment	84357129PMul2W
Angled Multiple Abutment	84357129PMulAngRA
Ti-BASE CAD-CAM Abutment	84357129CCPTiBRG
Ball-attached	84357129BallAttaBP
Provisional Abutment Ti PM	84357129PProvTiXS
Healing abutment	84357129TaCic3P
Abutment screw	84357129TorProt8Z
Micro Screw	84357129MicroTorS8

You can access the European Database on Medical Devices (EUDAMED) at the following link:

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